

Contec Products from HCT



peridox RTU
ANTIMICROBIAL SOLUTION
MOLD KILLING
Bactericidal, Virucidal,
Fungicidal, Sporicidal

Active Ingredients:
Sodium Hypochlorite 4.40%
Hydroxyethylamine 0.23%
Clear Ingredients 95.37%
TOTAL 100.00%

FIRST AID:

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor immediately for treatment advice. If swallowed, Call a poison control center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person. Have the product container or label with you when calling a poison control center or doctor for treatment.

IF SWALLOWED: Call a poison control center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person. Have the product container or label with you when calling a poison control center or doctor for treatment.

See back panel for additional precautionary statements.

EPA Reg. No. 6383-13-C-2; 67811-MQ-1; 71670-SC-1; 71670-SC-2;
EPA Est. No. 67720-MQ-1; 63838-AR-01; 67057-RC-1;
67057-RC-1; 63838-CA-01; 63838-AR-01; 67057-RC-1

1 U.S. Quart (0.9 liters)

PeridoxRTU is a trademark of Contec, Inc., Spartanburg, SC USA Patents Pending 2679-07/21

**KEEP OUT OF REACH OF CHILDREN
WARNING ADVERTENCIA**

If you do not understand the label, find someone to explain it to you in detail. Si usted no entiende la etiqueta, busque a alguien para que se la explique a usted en detalle. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person. Have the product container or label with you when calling a poison control center or doctor for treatment.

Emergency Telephone
Toll free 1-800-222-1222

















Proudly Providing Pharmacy
Solutions for Over 30 Years

Contec PreSaturated Wipes



Contec Wipes Designed that Exceed USP Standards
















Order No.	Description	Size	Qty
 HCPS6044	Contec Critical Site® Sterile wipes	4"x4"	1440 wipes/CS (60/pouch)
 HCPS2328IR	PROSAT® Sterile presaturated meltblown polypropylene wipes with 70% IPA/30% DI water	9"x11"	810 wipes/CS (30/pouch)
 HCPS0002IR	PROSAT® Sterile presaturated cellulose/polyester wipes with 70% IPA/30% DI water	9"x11"	810 wipes/CS (30/pouch)
 HCPS2328	PROSAT® presaturated meltblown polypropylene wipes with 70% IPA/30% DI water	9"x11"	810 wipes/CS (30/pouch)
 HCPS0002	PROSAT® presaturated cellulose/polyester wipes with 70% IPA/30% DI water	9"x11"	810 wipes/CS (30/pouch)
 PS-911EB	PROSAT® Sterile presaturated meltblown polypropylene wipes with 70% IPA/30% DI water	9"x11"	1440 wipes/CS (30/pouch)
 PS-911	PROSAT® presaturated meltblown polypropylene wipes with 70% IPA/30% DI water	9"x11"	1440 wipes/CS (30/pouch)
 PS-911-12	PROSAT® presaturated meltblown polypropylene wipes with 70% IPA/30% DI water	9"x11"	360 wipes/CS (30/pouch)
 PSPS0076	PROSAT® Sterile Polynit Heatseal Half-Fold Wipes with 70% IPA/30% DI water	9"x9"	800 wipes/CS (20/pouch)
 HCTB3315IR	Contec Healthcare TB1-3300™ disinfectant Sterile wipes	9"x11"	810 wipes/CS (30/pouch)
 HCTB3315	Contec Healthcare TB1-3300™ disinfectant wipes	9"x11"	810 wipes/CS (30/pouch)
 PSC20005	PROSAT® Sterile Theta™ wipes, lightweight nonwoven cellulose/polyester presaturated 70% IPA/30% DI water	9"x11"	1400 wipes/CS (28/pouch)
 21221	Contec PREempt RTU disinfectant, meltblown polypropylene wipes with Accelerated Hydrogen Peroxide	6"x7"	12/CS (160/canister)
 2B221	PREempt Plus Meltdown Wipes, polypropylene wipes with Accelerated Hydrogen Peroxide	6"x7"	12/CS (160/canister)









Contec Low-Lint, Sterile/Non-Sterile Dry Wipes



Order No.	Description	Size	Qty
 HCC20001	Absorbent, disposable, low-linting Hand Drying wipes.	17"x9"	900 wipes/CS (60/pouch)
 9451	Wall mounted stainless steel dispenser for HCC20001 hand drying wipes	N/A	1/CS
 NWPZ0001	Amplitude™ Kappa™ Sterile LE wipes lyocell/polyester blend nonwoven fabric	9"x9"	1000 wipes/CS (30/pouch)
 NWPZ0002	Amplitude™ Kappa™ Sterile LE wipes lyocell/polyester blend nonwoven fabric	12"x12"	1000 wipes/CS (30/pouch)
 AMEP0001	Amplitude™ Epsilon™ wipes cellulose/polyester nonwoven wipes	9"x9"	3600 wipes/CS (30/pouch)
 AMEP0002	Amplitude™ Epsilon™ wipes cellulose/polyester nonwoven wipes	12"x12"	3600 wipes/CS (30/pouch)
 AMEC0002	Amplitude™ EcoCloth™ wipes cellulose/polyester nonwoven wipes	6"x6"	6000 wipes/CS (30/pouch)
 AMEC0003	Amplitude™ EcoCloth™ wipes cellulose/polyester nonwoven wipes	9"x9"	3600 wipes/CS (30/pouch)
 AMEC0004	Amplitude™ EcoCloth™ wipes cellulose/polyester nonwoven wipes	12"x12"	2100 wipes/CS (30/pouch)
 AMEC0005	Amplitude™ EcoCloth™ wipes cellulose/polyester nonwoven wipes	18"x18"	750 wipes/CS (10/pouch)
 AMSI0001	Amplitude™ Sigma™ Wipes Cellulose/polyester hydroentangled nonwoven wipes	9"x9"	3600 wipes/CS (30/pouch)
 C2-99IR/25	SterileSorb™ wipes Sterile hydroentangled cellulose/polyester wipes	9"x9"	1050 wipes/CS (30/pouch)
 C2-1212IR	SterileSorb™ wipes Sterile hydroentangled cellulose/polyester wipes	12"x12"	1050 wipes/CS (30/pouch)

Contec Dry Wipes

Order No.	Description	Size	Qty
 LWPS0006	Polynit Heatseal Sterile wipes textured polyester interlock knit	9"x9"	1200 wipes/CS (25/pouch)
 LWMM0300	Microfiber PET/PP wipes Cleanroom laundered, bulk	10"x18"	400 wipes/CS (100/pouch)
 LWMM0400	Microfiber PET/PP wipes Cleanroom laundered, bulk	10"x18"	1 kit/CS
 HCPW3001	Microfiber PET/PP wipes Cleanroom laundered, bulk	10"x12"	700 wipes/CS (25/pouch)
 HCPW3003	Microfiber PET/PP wipes Cleanroom laundered, bulk	10"x18"	250 wipes/CS (25/pouch)
 85337	Sanotex® environmental surface wipes Half fold	10"x12"	1800 wipes/CS (60/roll)
 85338	Sanotex® environmental surface wipes Half fold	10"x12"	320 wipes/CS (40/roll)
 85339	Sanotex® environmental surface wipes	6"x8"	900 wipes/CS (150/roll)
 85340	Sanotex® environmental surface wipes Half fold	10"x12"	700 wipes/CS (10/roll)
 85341	Sanotex® environmental surface canister Empty self-closing dispenser, reusable	N/A	6 canisters/CS



Contec Cleaning and Disinfection Chemicals







Specially-formulated sporicidal agents, disinfectants, and chemicals for cleaning and disinfecting sterile/critical pharmacy environments; including 70% Isopropanol, hydrogen peroxide-based solutions, to fast-acting sporicidal disinfectants and broad-spectrum daily disinfectants.



Order No.	Description	Size	Qty
 HC85335	PeridoxRTU broad spectrum sporicidal disinfectant EPA Reg. No. 8383-13 Flip-top cap	32 oz.	6/CS
 HC85335IR	PeridoxRTU, Sterile broad spectrum sporicidal disinfectant EPA Reg. No. 8383-13, Flip-top cap	32 oz.	6/CS
 HC85335IR	PeridoxRTU broad spectrum sporicidal disinfectant EPA Reg. No. 8383-13	1 gal.	6/CS
 HC85336IR	PeridoxRTU, Sterile broad spectrum sporicidal disinfectant EPA Reg. No. 8383-13	1 gal.	6/CS
 HCTB3315IR	Contec Healthcare TB1-3300 disinfectant Sterile wipes	9"x11"	810 wipes/CS (30/pouch)
 HCTB3315	Contec Healthcare TB1-3300 disinfectant wipes	9"x11"	810 wipes/CS (30/pouch)
 HCTB3300	Contec Healthcare TB1-3300 disinfectant Single bagged, push-pull cap	32 oz.	12/CS
 HCTB3300IR	Contec Healthcare TB1-3300 disinfectant Sterile , double bagged, push-pull cap	32 oz.	12/CS
 HCTB3305	Contec Healthcare TB1-3300 disinfectant Single bagged	1 gal.	4/CS
 HCTB3305IR	Contec Healthcare TB1-3300 disinfectant Sterile , double bagged	1 gal.	4/CS
 SB167030IR	Contec Sterile 70% isopropanol Trigger sprayer	16 oz.	12/CS
 SB327030IR	Contec Sterile 70% isopropanol Trigger sprayer	32 oz.	12/CS
 SB1287030IR	Contec Sterile 70% isopropanol	1 gal.	4/CS

Contec Critical Pharmacy Products

Contec Cleaning and Disinfection Chemicals (cont.)



Order No.	Description	Size	Qty
 HCFT7030IR	Contec Sterile 70% isopropanol Flip-top cap	16 oz.	12/CS
 HCFT7030IR-32	Contec Sterile 70% isopropanol Flip-top cap	32 oz.	12/CS
 21101	Contec PREempt RTU Disinfectant With Accelerated Hydrogen Peroxide, capped	32 oz.	12/CS
 21105	Contec PREempt RTU Disinfectant With Accelerated Hydrogen Peroxide, capped	1 Gallon	4/CS
 2B101	Contec PREempt Plus One-Step Disinfectant Cleaner	32 oz.	12/CS
 2B105	Contec PREempt Plus One-Step Disinfectant Cleaner	1 Gallon	4/CS

Contec Low-Lint Mops for Cleanrooms & More







Systems and accessories for low-lint mopping that meet and exceed USP Standards, including mop frames, handles, wall washing systems, single-use disposable microfiber mop pads, and more

















Contec Cleaning Pads / Mop Covers (One-Time Use, Disposable)

Order No.	Description	Size	Qty
 PSME0001	EasyReach™ cleaning pads Quilted polyester presaturated with 70% IPA	2.75" x 7.75"	64/CS
HCMT0011	EasyReach™ cleaning pads with ultrasonic pockets. Quilted polyester	2.75" x 7.75"	60/CS
HCMT0015	EasyReach™ cleaning pads with ultrasonic pockets. Quilted polyester, Sterile	2.75" x 7.75"	60/CS
 HCPM5001	MicroCinch™ Mop Microfiber mop cover with pockets. Standard weight	5" x 15.5"	120/CS
HCPM5051	MicroCinch™ Mop Microfiber mop cover with pockets. Standard weight	5" x 15.5"	70/CS

Contec Mop Heads (One-Time Use, Disposable)

Order No.	Description	Size	Qty
 HCKM3001	Klean Max™ mop head, small	8" x 4" x 1"	48/CS
HCKM3002	Klean Max™ mop head, medium	13" x 4" x 1"	28/CS
HCKM3003	Klean Max™ mop head, large	15" x 4" x 1"	24/CS
HCKM3051	Klean Max™ mop head, small. Sterile	8" x 4" x 1"	48/CS
HCKM3052	Klean Max™ mop head, medium. Sterile	13" x 4" x 1"	28/CS
HCKM3053	Klean Max™ mop head, large. Sterile	15" x 4" x 1"	24/CS
 HCKM3005	Klean Max™ mop head, small. Sealed edge	9.5" x 5" x 0.375"	40/CS
HCKM3004	Klean Max™ mop head, medium. Sealed edge	15" x 5" x 0.375"	16/CS
HCKM3055	Klean Max™ mop head, small Sealed edge, Sterile	9.5" x 5" x 0.375"	40/CS
HCKM3054	Klean Max™ mop head, medium Sealed edge, Sterile	15" x 5" x 0.375"	24/CS
 EZMP0350	EasyCurve™ mop head. Microfiber individually packaged, Sterile	14" x 4.3" x .75"	1/CS
EZMP0300	EasyCurve™ mop head. Microfiber	14" x 4.3" x .75"	1/CS
 TASK0100	QuickClean foam mop head Laminated with polyester fabric	N/A	32/CS
 HCPM2002	Premira II White Microfiber Pads	N/A	200/CS
 BLVKMOPIR	EasySat™ Mop Head, Sterile Knitted polyester/foam mop	21.5" x 18.25" x 4.75"	28/CS

Contec Mops & Handles

Contec Mop Handles / Mop Frames / Tools				
Order No.	Description	Size	Qty	
	2642 QuickConnect™ anodized aluminum handle (Compatible with EasyReach™)	16" -30"	1/CS	
	2643 QuickConnect™ anodized aluminum handle (Compatible with EasyReach™)	50" -92"	1/CS	
	2644 QuickConnect™ stainless steel handle (Compatible with EasyReach™ & EasyCurve™)	16" -30"	1/CS	
	2645 QuickConnect™ stainless steel handle (Compatible with EasyCurve™)	50" -92"	1/CS	
	2725E QuickConnect™ stainless steel handle (Compatible with MicroCinch™, EasyCurve™, Klean Max™, & QuickClean™)	60"	1/CS	
	2746S QuickConnect™ handle Fiberglass with stainless steel tip	40" -70"	1/CS	
	2685LF QuickConnect™ EasyReach™ mop frame, layflat. Stainless steel	7.5"	1/CS	
	2665LF EasyReach™ Cleaning Tool, layflat, lightweight Stainless steel	16.5"	1/CS	
	HCPH8007 MicroCinch™ flat hinged mop frame Stainless steel	16"x4.25"	1/CS	
	2651SS QuickConnect™ mop frame, small For Klean Max™ small mop head	7.5"	1/CS	
	2652LF QuickConnect™ mop frame, layflat, medium For Klean Max™ medium and large mop heads	11.5"	1/CS	
	2740 MicroCinch™ flat hinged mop frame Blue polypropylene	16"x4.25"	1/CS	
	2745 QuickClean flat mop head frame Blue polypropylene	16"	1/CS	
	2724 EasyCurve™ pivoting mop frame Stainless steel	2.75"x5.25"x14.5"	1/CS	



Cleanroom garb for critical and USP-compliant environments, including gowns, frocks, garb, sterile personal protective equipment (PPE), and more



Contec CritiGear™ Sterile Nitrile Gloves (One-Time Use, Disposable)

- Contec CritiGear Sterile Nitrile Gloves have been specifically designed for pharmacy environments regulated by USP <797> and USP <800>.
- Sterility Assurance Level 10^{-6} | Tested per ASTM D6978-05

Order No.	Description	Packaging	Qty
HCGS0002	X-Small	1 pair/pack	200 pair/CS
HCGS0012	Small	1 pair/pack	200 pair/CS
HCGS0022	Medium	1 pair/pack	200 pair/CS
HCGS0032	Large	1 pair/pack	200 pair/CS
HCGS0042	X-Large	1 pair/pack	200 pair/CS
HCGS0052	2X-Large	1 pair/pack	200 pair/CS

Contec Shoe Covers (One-Time Use, Disposable)

Order No.	Description	Size	Qty
 SCR645L-150	SCR645 series shoe covers. Blue	Universal	300 each/CS
SCR645XL-150	SCR645 series shoe covers. Blue	X-Large	300 each/CS
SCR645XXL-150	SCR645 series shoe covers. Blue	2X-Large	300 each/CS
 SCPL500	SCPL series compressed polyethylene shoe covers. White	Universal	1000 each/CS
 LKBT0001-S/M-IR	CleanMax™ STERILE 19" Boot Cover Non-Skid PVC vinyl sole, tunneled elastic top	Small/Medium	50/CS

Contec Garb & PPE (cont.)



Bouffant Cap 24" - White

- Strong, Durable, & Lightweight with No Latex Proteins
- Contamination Control
- Individually Inspected for Stitching, Sizing, Material and Workmanship

Order No.	Description	Qty
NBC2400	24"	1000/CS



Beard Cover White Spunbonded Poly - White

- Strong, Durable, & Lightweight with No Latex Proteins
- Contamination Control
- Individually Inspected for Stitching, Sizing, Material and Workmanship

Order No.	Description	Qty
HBC1000	Approx. Length: 38.1cm (15 in.)	1000/CS



Bouffant Cap 24" - White

- Strong, Durable, & Lightweight with No Latex Proteins
- Contamination Control
- Individually Inspected for Stitching, Sizing, Material and Workmanship

Order No.	Description	Qty
GOGLIR	24"	1000/CS



Contec Cleanroom Mask - White

- Earloop Face Mask
- Strong, Durable, & Lightweight
- Cleanroom Rated

Order No.	Description	Qty
MAS-9010	White	500/CS



Contec CritiGear™ Cleanroom Frock with Thumb Loop

Order No.	Size	Qty
HCGA0002	X-Small	30/CS
HCGA0012	Small	30/CS
HCGA0022	Medium	30/CS
HCGA0032	Large	30/CS
HCGA0042	X-Large	30/CS
HCGA0052	2X-Large	30/CS
HCGA0062	3X-Large	30/CS
HCGA0072	4X-Large	30/CS
HCGA0082	5X-Large	30/CS



Contec CritiGear™ Chemo Hazardous Drug Gowns with Thumb Loop

Order No.	Size	Qty
HCGA0800-S/M	Small/Medium	30/CS
HCGA0800-L/XL	Large/X-Large	30/CS
HCGA0800-2XL/3XL	2X-Large/3X-Large	30/CS
HCGA0800-4XL/5XL	4X-Large/5X-Large	30/CS
HCGA0800-6XL/7XL	6X-Large/7X-Large	30/CS



CleanMax™ Coveralls with Hood, Tunneled Wrists, & Ankles

Order No.	Size	Qty
CVHD0022	Medium	25/CS
CVHD0032	Large	25/CS
CVHD0042	X-Large	25/CS
CVHD0052	2X-Large	25/CS
CVHD0062	3X-Large	25/CS

Contec Garb & PPE (cont.)

STERILE



CleanMax™ garments feature bound seams, which are sewn with an additional outer binding to increase seam strength and provide better barrier protection.



CleanMax™STERILE Coveralls with Hood, Tunneled Wrists, & Ankles

Order No.	Size	Qty
CVHD0022IR	Medium	25/CS
CVHD0032IR	Large	25/CS
CVHD0042IR	X-Large	25/CS
CVHD0052IR	2X-Large	25/CS
CVHD0062IR	3X-Large	25/CS



CVRL

CleanMax™ Coveralls with Tunneled Wrists & Ankles

Order No.	Size	Qty
CVRL0022	Medium	25/CS
CVRL0032	Large	25/CS
CVRL0042	X-Large	25/CS
CVRL0052	2X-Large	25/CS
CVRL0062	3X-Large	25/CS

STERILE



CVRL

CleanMax™STERILE Coveralls with Tunneled Wrists & Ankles

Order No.	Size	Qty
CVRL0022IR	Medium	25/CS
CVRL0032IR	Large	25/CS
CVRL0042IR	X-Large	25/CS
CVRL0052IR	2X-Large	25/CS
CVRL0062IR	3X-Large	25/CS



Contec Critical Environment Supplies


Critical supplies for your critical environment, including prep and tack mats, disposable mesh foam sponges, and an essential hazardous drug spill kit



Clean Scrub™ Sponge (Disposable mesh sponge with foam core)

- Textured polyester mesh removes heavily soiled contamination from surfaces
- Compatible with a wide range of disinfectants

Order No.	Size	Qty
VKSC0004	5"x3"x0.6"	48 each/CS



Contec Healthcare Essential HD Spill Toolkit (for up to 1000mL)

- Activated charcoal spill pads begin chemical deactivation of spills, absorbs, and renders liquid drug unrecoverable, enhancing the safety of workers absorbing up to 1000mL per pad
- Integrated Chemo Bio-Wipe bag safely contains macro spill components for later disposal in RCRA* container
- Plastic scoop and rigid container conveniently facilitate powder spill cleanup and containment

Order No.	Qty
HCSK0001	2 kits/CS



Irradiated Bucket Liners

- 4 mil polyethylene; fits 5 gallon, 6.5 gallon, 7 gallon
- 2/bag, 20 bags/case


Order No.	Size	Qty
BKTBAG	34"x24"x.004"	40/CS



Contec Prep Mat - Contaminaton Control Mat

- Low-lint construction; absorbent and barrier layers
- Lot number on inner bags for traceability
- Available irradiated and nonsterile

Order No.	Description	Size	Qty
HCPD0001	Contec prep mat	16"x21"	600 each/CS
HCPD0005	Contec prep mat. Irradiated	16"x21"	600 each/CS



Contec Tack Mat - Contaminaton Control Flooring

- Placed immediately before the entrance to controlled areas to reduce the introduction of contamination
- Thirty peelable layers for convenience
- Multiple sizes are available
- Frames are an optional way to avoid exposing the floor to adhesives

Order No.	Description	Size	Qty
TMF2436	Tacky mat frame	24"x36"	1 each/CS
TMF1836	Tacky mat frame	18"x36"	1 each/CS
TM0001	Tacky mat. Blue	24"x36"	120/CS
TM0006	Tacky mat. Blue	18"x36"	120/CS

Talk to Your Sales Rep About Product Programs & Availability