BiTrac® NIV

Non-Invasive Ventilation Masks

The BiTrac® NIV product family features a variety of interfaces that are designed to increase patient comfort and compliance in respiratory distress scenarios. The full-face and nasal masks are highly adjustable, and their cushion designs minimize pressure exerted on the bridge of the nose. Our total perimeter mask, the MaxShield™, avoids the nose and disperses pressure over a larger surface area while reducing claustrophobia for the patient. Our range of interchangeable elbows greatly reduces waste since you can easily remove and replace the mask's elbow when transitioning patients between standard ventilators and bi-level machines.

Interchangeable Elbows

Soft Silicone Cushion

Variety of sizes





Custom fit for Each Patient



Interchangeable Elbows



BiTrac™ Circuit

BiTrac® and MaxShield™ NIV Mask

Davit No.	Made	Ci=o	Interchangeable/Fired	Elhaur	Casa Oty
Part No.	Mask	Size	Interchangeable/Fixed	Elbow	Case Qty.
313-9013W	BiTrac Nasal	Adult Small	Fixed	Standard	10
313-9014W	BiTrac Nasal	Adult Medium	Fixed	Standard	10
313-9015W	BiTrac Nasal	Adult Large	Fixed	Standard	10
313-9016W	BiTrac Nasal	Adult Extra Large	Fixed	Standard	10
313-9001W	BiTrac Full Face	Adult Small	Fixed	Anti-Asphyxia	10
313-9002W	BiTrac Full Face	Adult Medium	Fixed	Anti-Asphyxia	10
313-9003W	BiTrac Full Face	Adult Large	Fixed	Anti-Asphyxia	10
313-9004W	BiTrac Full Face	Adult Extra Large	Fixed	Anti-Asphyxia	10
313-9009W	BiTrac Full Face	Adult Small	Fixed	Standard	10
313-9010W	BiTrac Full Face	Adult Medium	Fixed	Standard	10
313-9011W	BiTrac Full Face	Adult Large	Fixed	Standard	10
313-9012W	BiTrac Full Face	Adult Extra Large	Fixed	Standard	10
313-9501	BiTrac Select Full Face	Adult Small	Interchangeable	Anti-Asphyxia	10
313-9502	BiTrac Select Full Face	Adult Medium	Interchangeable	Anti-Asphyxia	10
313-9503	BiTrac Select Full Face	Adult Large	Interchangeable	Anti-Asphyxia	10
313-9504	BiTrac Select Full Face	Adult Extra Large	Interchangeable	Anti-Asphyxia	10
313-9506	BiTrac Select Full Face	Adult Small	Interchangeable	Standard	10
313-9507	BiTrac Select Full Face	Adult Medium	Interchangeable	Standard	10
313-9508	BiTrac Select Full Face	Adult Large	Interchangeable	Standard	10
313-9509	BiTrac Select Full Face	Adult Extra Large	Interchangeable	Standard	10
313-9555	BiTrac MaxShield Select	Adult Extra Large	Interchangeable	Anti-Asphyxia with Leak	10
313-9552	BiTrac MaxShield Select	Adult Large	Interchangeable	Anti-Asphyxia with Leak	10
313-9554	BiTrac MaxShield Select	Adult Small	Interchangeable	Anti-Asphyxia with Leak	10
313-9556	BiTrac MaxShield Select	Adult Extra Large	Interchangeable	Standard	10
313-9553	BiTrac MaxShield Select	Adult Large	Interchangeable	Standard	10
313-9551	BiTrac MaxShield Select	Adult Small	Interchangeable	Standard	10
313-9057W	BiTrac MaxShield	Adult Extra Large	Fixed	Anti-Asphyxia	10
313-9052W	BiTrac MaxShield	Adult Large	Fixed	Anti-Asphyxia	10
313-9054W	BiTrac MaxShield	Adult Small	Fixed	Anti-Asphyxia	10
313-9056W	BiTrac MaxShield	Adult Extra Large	Fixed	Standard	10
313-9053W	BiTrac MaxShield	Adult Large	Fixed	Standard	10
313-9055W	BiTrac MaxShield	Adult Small	Fixed	Standard	10
1.1.7000				2.2	

BiTrac® Accessories

Part No.	Description	Case Qty.
313-9593	Interchangeable Elbow - Anti-Asphyxia with Leak Elbow	10
313-9594	Interchangeable Elbow - Anti-Asphyxia Nebulizer Elbow	10
313-9595	Interchangeable Elbow - Standard Nebulizer Elbow	10
313-9596	Interchangeable Elbow - Anti-Asphyxia Bronchoscopy Elbow	10
313-9597	Interchangeable Elbow - Standard Bronchoscopy Elbow	10
313-9598	Interchangeable Elbow - Anti-Asphyxia Elbow	10
313-9599	Interchangeable Elbow - Standard Elbow	10
313-9074	BiTrac NIV Circuit (Bilevel/CPAP Application), 72" Smooth ID Tubing, swivel connector with exhalation port, Cap Plug, 22mmF x 22mmF Adapter	10
313-9075	BiTrac NIV Circuit (Bilevel/CPAP Application), 72" Smooth ID Tubing with attached filter, swivel connector with exhalation port, Cap Plug, 22mmF x 22mmF Adapter	10
313-9078	ValueTrac NIV Circuit (Bilevel/CPAP Application), 72" Corrugated Tubing and swivel connector with exhalation port and optional 22mmF x 22mmF Adapter	10