



Tamper Evident Products from HCT



Proudly Providing Pharmacy
Solutions for Over 30 years



HCT Tamper Evident Products Increase Security and Value!

Tamper evident caps and port clamps from HCT ensure product integrity from the moment the drug leaves the PEC until it is administered by clinicians. Show that your facility takes security seriously as these products increase your accreditation and provide significant value.

Hospitals and facilities bear the cost of diverted drugs, and the government can fine up to **\$4M** for **inadequate safeguards**. Reduce the risk of being compromised either accidentally or intentionally.

Custom Colors Available with Volume Agreements

Talk to your sales rep about product programs and availability

Overall Pros of HCT Tamper Evident Syringe Caps & Port Clamps:

- Active, low-cost deterrent to diversion
- Easy to use
- Provides a physical barrier to the opportunity for misuse while enhancing medication safety
- Different color caps allow the pharmacy to color-code by drug classes, processes, or however they choose
- Reduces risk of contamination, infection, and incorrect dosing for patient
- Ideal for the majority of hospitals continuing to implement drug awareness initiatives
- Reduces waste & associated cost
- Simple to apply
- Streamlines processes while supporting consistent patient care
- Not required by the DEA, but strengthens and supplements DEA and USP guidelines
- Increases patient safety by preventing leakage and ensuring they receive proper dosage
- Strengthens <797> compliance





The Benefits of Tamper Evident Caps

- Simply twists onto the syringe
- Syringes with intact caps may be returned to the pharmacy
- Maintains aseptic technique and saves time when compounding several syringes



Color Inside Only



Color Inside/Outside

Custom Colors Available with Volume Agreements

Order No.	Color
-----------	-------

Order No.	Color
-----------	-------

Caps for Small Volume Compounders

Easy open packaging allows for access to 2 caps at a time, reducing risk of contamination to other caps

10 per Tray | 500/Case

HCT10R	Red
HCT10W	White
HCT10B	Blue
HCT10K	Black

HCT10RR	Red/Red
HCT10WW	White/White
HCT10BB	Blue/Blue
HCT10KK	Black/Black



Caps for High Volume Manual Compounders

Increase throughput and reduce sanitization time and materials

100 Pack (10 per Tray) | 1000/Case

HCT100R-5B	Red
HCT100W-5B	White
HCT100B-5B	Blue
HCT100K-5B	Black

HCT100RR-5B	Red/Red
HCT100WW-5B	White/White
HCT100BB-5B	Blue/Blue
HCT100KK-5B	Black/Black

All trays placed in sterile overwrap

Caps for Automated Syringe Fillers

Easy to open notch to work quickly with automated syringe filling machines

Bulk Sterile Bags 500 per Bag | 2 Bags per CS

HCT500R-5B	Red
HCT500W-5B	White
HCT500B-5B	Blue
HCT500K-5B	Black

HCT500RR-5B	Red/Red
HCT500WW-5B	White/White
HCT500BB-5B	Blue/Blue
HCT500KK-5B	Black/Black



The Benefits of Tamper Evident Port Clamps

- An active, low-cost deterrent to diversion
- Unadministered bags can be returned to the pharmacy reducing shortages
- Increases IV bag security while reducing waste



Custom Colors Available with Volume Agreements

Order No.	Description	Color	Qty
-----------	-------------	-------	-----


Tamper Port Clamps and Their Compatibility by Brand

Easy to install for guaranteed security!

HCT50-BXR	Tamper Evident Additive Port Clamp for IV Bags Compatible with Baxter Viaflex, and Intravia Bags	Red	50/CS
HCT50-MXR	Tamper Evident Additive Port Clamp for IV Bags Compatible with Metrix, and HCT Bags	Red	50/CS
HCT50-HOSXR	Tamper Evident Additive Port Clamp for IV Bags Compatible with Hospira Bags	Red	50/CS
HCT50-BBXR	Tamper Evident Additive Port Clamp for IV Bags Compatible with B.Braun Excel Bags	Red	50/CS
HCT50-BBPABR	Tamper Evident Additive Port Clamp for IV Bags Compatible with B.Braun Pab Bags	Red	50/CS



HCT50-HOSXR Red Shown



HCT50-MXR Red Shown