

HMEs and Filters Available in Sterile or Clean Packaging

Ask your local HCT rep about how to secure your inventory with our stocking agreements for HMEs, HMEs and Filters, HEPA Filters, and Electrostatic Filters.

Bacterial Viral Filters & More!



Product Codes Clean:
313-435

TV Range (ml):

Weight (grams):

Resistance (cm H₂O):

Dead Space (ml):

Box/Case Qty.:

Bacterial Efficiency %: >99.999
Viral Efficiency %: >99.99

>125

22.5

1.0,60 LPM

30.0

20/200

Brands You Trust at Your Fingertips, Ask Your HCT Rep for Details Today!

