

## Passive Bilevel 22mm Breathing Systems

Intersurgical have a complete range of Flextube™ bilevel breathing systems compatible with turbine driven bilevel machines that require a single limb breathing system.

Our bilevel products offer breathing system solutions with CO<sub>2</sub> leak port for use with standard NIV (Non-invasive Ventilation) masks, and systems without CO<sub>2</sub> leak port for use with ported NIV masks. Breathing systems without proximal airway pressure monitoring lines are also available for those machines that do not require this feature. Breathing systems in this range are suitable for adult and paediatric patients where specified.

- 22mm adult/paediatric
- Smoothbore bilevel breathing system with swivel patient connector
- CO<sub>2</sub> leak port and monitoring line
- ≥1.6m



| Order No. | Description                            | Qty   |
|-----------|--|-------|
| 5805000   | Passive Bilevel 22mm Breathing Circuit | 35/CS |

