





Disposable, Single-Use HCT Dial-A-Flow

Talk to your sales rep about product programs and availability.
 Orders will ship the same day if placed before 4 PM EST.
 All Dial-A-Flow models feature 28.2 ml priming.



Order No.	Features	Qty
HCT100DFSDM	 <ol style="list-style-type: none"> 1. Universal ABS+PVC Spike & Drip Chamber 2. ABS Upper Pinch Clamp 3. Dial-Type ABS Flow Regulator 4. Rotating Make Luer Lock 5. Non-DEHP PVC Tubing 	50/CS
HCT100DFSDM-F	 <ol style="list-style-type: none"> 1. Universal ABS+PVC Spike & Drip Chamber 2. ABS Upper Pinch Clamp 3. Dial-Type ABS Flow Regulator 4. Rotating Make Luer Lock 5. Non-DEHP PVC Tubing 6. Acrylic 0.2µ Filter 	50/CS
HCT100DFSDM-FNY	 <ol style="list-style-type: none"> 1. Universal ABS+PVC Spike & Drip Chamber 2. ABS Upper Pinch Clamp 3. Dial-Type ABS Flow Regulator 4. Rotating Make Luer Lock 5. Non-DEHP PVC Tubing 6. Acrylic 0.2µ Filter 7. Needle-Free Y-Site 	50/CS
HCT100DFSDM-NY	 <ol style="list-style-type: none"> 1. Universal ABS+PVC Spike & Drip Chamber 2. ABS Upper Pinch Clamp 3. Dial-Type ABS Flow Regulator 4. Rotating Make Luer Lock 5. Non-DEHP PVC Tubing 6. Needle-Free Y-Site 	50/CS