



Pharmacy Products
from HCT



Proudly Providing Pharmacy
Solutions for Over 30 Years



Proudly Providing Pharmacy Solutions for Over 30 Years

HCT Products Promote Efficiency and Aseptic Technique

HCT is pleased to have partnered with leading pharmacy companies that align with our value based model. Please contact us to learn more about how our portfolio is delivering exceptional clinical outcomes and providing cost savings.

Talk to your sales rep about product programs and availability.

Table of Contents

Page 3-5	HCT Tamper Evident Syringe Caps & IV Bag Port Clamps
Page 6	HCT Admixture Products
Page 7-9	IV Bags
Page 10-11	HCT Peristaltic Pump Sets & Single-Use Systems (SUS)
Page 12-13	HCT Syringe Trays
Page 14-16	Filters
Page 17-19	HCT Personal Protective Equipment (PPE)
Page 20-24	QI Medical Quality Assurance Products
Page 25-37	Contec Products
Page 38-41	Medipak UVLI® Covers and Bags





HCT Tamper Evident Products Increase Security and Value!

Tamper evident caps and port clamps from HCT ensure product integrity from the moment the drug leaves the PEC until it is administered by clinicians. Show that your facility takes security seriously as these products increase your accreditation and provide significant value.

Hospitals and facilities bear the cost of diverted drugs, and the government can fine up to **\$4M** for **inadequate safeguards**. Reduce the risk of being compromised either accidentally or intentionally.

Custom Colors Available with Volume Agreements

Talk to your sales rep about product programs and availability

Overall Pros of HCT Tamper Evident Syringe Caps & Port Clamps:

- Active, low-cost deterrent to diversion
- Easy to use
- Provides a physical barrier to the opportunity for misuse while enhancing medication safety
- Different color caps allow the pharmacy to color-code by drug classes, processes, or however they choose
- Reduces risk of contamination, infection, and incorrect dosing for patient
- Ideal for the majority of hospitals continuing to implement drug awareness initiatives
- Reduces waste & associated cost
- Simple to apply
- Streamlines processes while supporting consistent patient care
- Not required by the DEA, but strengthens and supplements DEA and USP guidelines
- Increases patient safety by preventing leakage and ensuring they receive proper dosage
- Strengthens <797> compliance





HCT Tamper Evident Syringe Caps

The Benefits of Tamper Evident Caps

- Simply twists onto the syringe
- Syringes with intact caps may be returned to the pharmacy
- Maintains aseptic technique and saves time when compounding several syringes



Color Inside Only



Color Inside/Outside

Custom Colors Available with Volume Agreements

Order No.	Color
-----------	-------

Order No.	Color
-----------	-------

Caps for Small Volume Compounders

10 per Tray | 500/Case

Easy open packaging allows for access to 2 caps at a time, reducing risk of contamination to other caps

HCT10R	Red
HCT10W	White
HCT10B	Blue
HCT10K	Black

HCT10RR	Red/Red
HCT10WW	White/White
HCT10BB	Blue/Blue
HCT10KK	Black/Black



Caps for High Volume Manual Compounders

100 Pack (10 per Tray) | 1000/Case

Increase throughput and reduce sanitization time and materials

HCT100R-5B	Red
HCT100W-5B	White
HCT100B-5B	Blue
HCT100K-5B	Black

HCT100RR-5B	Red/Red
HCT100WW-5B	White/White
HCT100BB-5B	Blue/Blue
HCT100KK-5B	Black/Black

All trays placed in sterile overwrap

Caps for Automated Syringe Fillers

Bulk Sterile Bags 500 per Bag | 2 Bags per Case

Easy to open notch to work quickly with automated syringe filling machines

HCT500R-5B	Red
HCT500W-5B	White
HCT500B-5B	Blue
HCT500K-5B	Black

HCT500RR-5B	Red/Red
HCT500WW-5B	White/White
HCT500BB-5B	Blue/Blue
HCT500KK-5B	Black/Black

The Benefits of Tamper Evident Port Clamps

- An active, low-cost deterrent to diversion
- Unadministered bags can be returned to the pharmacy reducing shortages
- Increases IV bag security while reducing waste



Custom Colors Available with Volume Agreements

Order No.	Description	Color	Qty
-----------	-------------	-------	-----


Tamper Port Clamps and Their Compatibility by Brand

Easy to install for guaranteed security!

HCT50-BXR	Tamper Evident Additive Port Clamp for IV Bags (Compatible with Baxter Viaflex, and Intravia Bags)	Red	50/CS
HCT50-MXB	Tamper Evident Additive Port Clamp for IV Bags (Compatible with Metrix, and HCT Bags)	Blue	50/CS
HCT50-MXR	Tamper Evident Additive Port Clamp for IV Bags (Compatible with Metrix, and HCT Bags)	Red	50/CS
HCT50-HOSXR	Tamper Evident Additive Port Clamp for IV Bags (Compatible with Hospira Bags)	Red	50/CS
HCT50-BBXR	Tamper Evident Additive Port Clamp for IV Bags (Compatible with B.Braun Excel Bags)	Red	50/CS
HCT50-BBPABR	Tamper Evident Additive Port Clamp for IV Bags (Compatible with B.Braun Pab Bags)	Red	50/CS



HCT50-HOSXR Red Shown











HCT50-MXR Red Shown

HCT Admixture Products

HCT Admixture Products Promote Efficiency and Aseptic Technique



	Order No.	Description	Qty
	HCT1100	Fluid Dispenser Luer H Connector Female Luer Lock/Female Luer Slip	50/CS
	HCT1200	Fluid Dispenser Luer Connector Bulk Pack: HCT1200-B 25 single sterile per outer bag, 4 bags per case 100/CS Talk to your sales rep about product programs and availability	50/CS
	HCT1300	Fluid Dispenser Luer Connector (Male/Male Connector) Bulk Pack: HCT1300-B 25 single sterile per outer bag, 4 bags per case 100/CS Talk to your sales rep about product programs and availability	50/CS
	HCT2000	Mini-Transfer Spike Bulk Pack: HCT2000-B 25 single sterile per outer bag, 4 bags per case 100/CS Talk to your sales rep about product programs and availability	50/CS
	HCT2001	Chemo Dispensing Pin Bulk Pack: HCT2001-B 25 single sterile per outer bag, 4 bags per case 100/CS Talk to your sales rep about product programs and availability	50/CS
	HCT3000	Bag Spike Luer Adapter with Two Way Valve Bulk Pack: HCT3000-B 25 single sterile per outer bag, 4 bags per case 100/CS Talk to your sales rep about product programs and availability	50/CS
	N0650	Non-Vented Vial Access Spike with NIS-6P (20mm Multi-Dose Access Device)	50/CS
	9902R	Red Male/Female Dead End Caps	200/CS
	N9902B	Blue Male/Female Dead End Caps	200/CS



HCT IV Bag Offerings

Durable, single-use IV bags offer superior performance and cost savings. Avoid the risk of punctures or leaks with the EVA (Ethylene Vinyl Acetate) bag options.

HCT offers **better pricing** and your sales rep can provide samples upon request.

Additional Options Available with Minimum Quantity Agreements

Talk to your sales rep about product programs and availability

IV Bag Highlights:

- Promotes cost savings
- Transparent for improved admixture inspection
- Sizes from 50 mL to 4000 mL
- Durable & flexible
- Practical and sturdy carrying handles
- Port for safe medication injection through synthetic septum cap
- Bag spike adapter for secure connection to infusion set



EVA IV Bags

EVA IV Bag Benefits:

- Made with EVA (ethylene vinyl acetate)
- 100% PVC-, DEHP-, and latex-free
- Increased flexibility and resistance reduces risk of punctures or leaks
- Luer lock adapter for easy filling and infusing



Order No.	Description	Qty
Secure 2-Port EVA (Ethylene Vinyl Acetate) Sterile Fluid Pathway Bags		
66042	IV Bags - 50mL	50/CS
66040	IV Bags - 100mL	50/CS
66045	IV Bags - 250mL	50/CS
66050	IV Bags - 500mL	50/CS
66055	IV Bags - 1000mL	50/CS
66043	IV Bags - 50mL (Bulk packed in a double poly bag without individual packaging)	50/CS
66041	IV Bags - 100mL (Bulk packed in a double poly bag without individual packaging)	50/CS
Secure EVA (Ethylene Vinyl Acetate) Sterile Gravity Bags		
66100	IV Bags 2 Leg Detachable - 1000mL	50/CS
66120	IV Bags 3 Leg Detachable - 1000mL	50/CS
Secure 3-Leg EVA (Ethylene Vinyl Acetate) Sterile Fluid Pathway Bags		
66405	IV Bags - Male Screw Connector - 250mL	50/CS
66410	IV Bags - Male Screw Connector - 500mL	50/CS
66415	IV Bags - Male Screw Connector - 1000mL	50/CS
66425	IV Bags - Male Screw Connector - 2000mL	50/CS

Order No.	Description	Qty
Secure 3-Leg EVA (Ethylene Vinyl Acetate) Sterile Fluid Pathway Bags (cont.)		
66430	IV Bags – Male Screw Connector – 3000mL	50/CS
66435	IV Bags – Male Screw Connector – 4000mL	50/CS
66205	IV Bags – Female Screw Connector – 250mL	50/CS
66210	IV Bags – Female Screw Connector – 500mL	50/CS
66215	IV Bags – Female Screw Connector – 1000mL	50/CS
66225	IV Bags – Female Screw Connector – 2000mL	50/CS
66230	IV Bags – Female Screw Connector – 3000mL	50/CS
66235	IV Bags – Female Screw Connector – 4000mL	50/CS

HCT 2-Port EVA (Ethylene Vinyl Acetate) Sterile


HCT2050	IV Bags – 50mL Sterile Bulk Pack (Triple Pkg, 20/Pkg, 7pkg/CS)	140/CS
HCT20100	IV Bags – 100mL Sterile Bulk Pack (Triple Pkg, 20/Pkg, 6pkg/CS)	120/CS
HCT10500	IV Bags – 500mL Sterile Bulk Pack (Triple Pkg, 10/Pkg, 7pkg/CS)	70/CS

Ask your local sales rep about HCT's sterile EVA bags that increase throughput and reduce sanitization time for high volume compounders.



HCT Peristaltic Pump Sets & Single-Use Systems (SUS)

Compatible with Top Brands (as listed below)



Order No.	Description	Qty
	HCT11 Fluid Transfer Peristaltic Tubing Set with Vented Spike / Luer Lock (Compatible with Baxa/Baxter Products)	10/CS
	HCT21 Fluid Transfer Peristaltic Tubing Set Luer Lock / Luer Lock (Compatible with Baxa/Baxter Products)	10/CS
	HCT31A Fluid Transfer Peristaltic Tubing Extension Set with ABS Luer Lock to ABS Luer Lock, Low Priming Volume (Compatible with Baxa/Baxter Products)	10/CS
	HCT86 Fluid Transfer Peristaltic Tubing Extension Set with Male Luer Lock to Female Luer Lock (Compatible with Baxa/Baxter Products)	50/CS
	HCT87A Fluid Transfer Peristaltic Tubing Repeater Extension Set with ABS Luer Lock to ABS Female Luer Lock (Compatible with Baxa/Baxter Products)	50/CS
	HCT29 Fluid Transfer Peristaltic Tubing Set Add-On with Vented Spike / Female Luer Lock (Compatible with Baxa/Baxter Products)	50/CS
	HCT550 Fluid Transfer Peristaltic Tubing Set Add-On with Male Luer Lock to Male Luer Lock - 60" (Compatible with Baxa/Baxter Products)	50/CS
	HCTPA01 PharmAssist Fluid Transfer Peristaltic Tubing Set - Single with Spike / Luer Lock, Low Priming Volume (Compatible with ICU PharmAssist)	10/CS
	HCTPA02 PharmAssist Fluid Transfer Peristaltic Tubing Set - Double with Spike / Luer Lock, Low Priming Volume (Compatible with ICU PharmAssist)	10/CS
	HCTPA03 PharmAssist Fluid Transfer Peristaltic Tubing Set - Triple with Spike / Luer Lock, Low Priming Volume (Compatible with ICU PharmAssist)	10/CS


HCT Can Meet Your Complex Configuration Needs

HCT is available for consultation to meet your design requirements.



Let Us Show You Our Customization Capabilities, Contact Your Local Sales Rep

Below is a custom peristaltic pump set & single-use system (SUS) that HCT created for a client.




No.	Park Description	Material	Remarks
1	3L EVA BAG	SANDED	30um
2	COEXTRUSION TUBE+MALE LUER	PVC(NON DEHP)&EVA+ABS	ID:6.5MM / OD:8.5MM / LENGTH:3CM LUER:5.0MM
3	FEMALE LUER (X2)	ABS	6.4MM
4	SILICONE TUBE (X2)	SILICONE	ID:6.4MM OD:9.5MM /LENGTH:6.4CM
5	CONNECTOR	PP	7.9MM x 6.4MM
6	SILICONE TUBE	SILICONE	ID:4.8MM OD:8.0MM /LENGTH:110CM
7	MALE LUER	PC	6.4MM
8	H CONNECTOR	ABS	BLUE
9	CONNECTOR	PP	7.9MM x 14.2MM
10	SILICONE TUBE	SILICONE	ID:12.7MM OD:19.05MM /LENGTH:10.16CM
11	CAPSULE FILTER1	707 (PP housing PES 0.22 Barbed hose ends)	Filtration area : 570cm ²
12	SILICONE TUBE (X2)	SILICONE	ID:6.4MM OD:11.2MM LENGTH:7CM
13	CAPSULE FILTER2	L707 (PP housing PES 0.22 Barbed hose ends)	Filtration area: 2041cm ²
14	FLANGE	PP	
15	SILICONE GASKET	SILICONE	
16	CLAMP	ABS	
17	END CAP	PP	
18	Metal Clamp (2X)	304	(0010.5-505R) 105
19	Metal Clamp (2X)	304	(0012.3-706R) 123
20-23	Metal Clamp (1X)	304	
24	Metal Clamp (3X)	304	(0015.7-706R) 157
25	Metal Clamp (1X)	304	(0014.2-706R) 142
26	PE BAG WITH TIE (2X)	PE+NYLONE	

Compatible with Brands Like:

- Bausch
- Cole-Parmer
- Flexicon
- Fisherbrand
- & More

HCT/BD Syringe Trays - Bulk Packaged

BD Syringe Information



Example of 5 Trays per Outer Bag

HCT's custom bulk packed BD Syringe trays are designed for high volume compounders. All syringes are laid unidirectionally in each tray, multiple trays are packaged in a poly overwrap and sterilized to ensure all contents inside the poly bag are sterile including the trays themselves. This packaging configuration allows for these syringes to move through different ISO classes in a highly-efficient way, reducing labor and sanitization materials.

- An additional AQL inspection is taken of the bulk non-sterile syringes including visual verification of any printing errors in the syringe graduations. It also includes visual inspection of contamination on the syringes, such as visible particles on the surfaces of the syringes.
- An additional wipe-down occurs upon receiving the bulk syringes in the bulk poly bags to remove particles and other materials.
- A daily visual inspection occurs dedicated to identifying foreign particles within the Cleanroom setting.

Packaging Summary

Syringe Material Information:

Barrels, Plunger Rods, Tip Shields: Polypropylene

Stopper: Synthetic Elastomer

Lubricant: Silicone

Syringe Tray Information:

All syringes oriented the same direction inside tray.

Tray Material: High Impact Polystyrene

Tray Lid Material: Coated Tyvek for breathable easy peel

Bag Material: 13" x 24" 4mil Polyethylene with 3" Tyvek header

Sterilization and Shelf Life:

Modality: Ethylene Oxide Sterilization following EN ISO11135-1. ETO residues are within applicable regulations


Shelf Life: 3 years

Bulk Packaged BD Syringes Ordering Details



Order No.	Description	Syringes per Tray	Trays per Outer Bag	Total Syringes per Case
HCTBD1ML-ST	1mL Slip Tip	50	5	1250/CS
HCTBD1ML-LL	1mL Luer Lock	50	5	1250/CS
HCTBD3ML-LL	3mL Luer Lock	50	5	1250/CS
HCTBD5ML-LL	5mL Luer Lock	30	5	750/CS
HCTBD10ML-LL	10mL Luer Lock	25	5	625/CS
HCTBD20ML-LL	20mL Luer Lock	TBD	TBD	TBD
HCTBD30ML-LL	30mL Luer Lock	TBD	TBD	TBD
HCTBD60ML-LL	60mL Luer Lock	TBD	TBD	TBD

Cardinal Single-Use Syringe Trays



- Sterile, single-use and bulk packaged syringes for pre-filling and batch filling procedures
- Not manufactured with DEHP or natural rubber latex

Order No.	Description	Tip Type	Syringes per Tray	Packs per Tray	Qty
8881501459	1mL Sterile Syringe Tray	Luer-Slip	25	40	1000/CS
8881513207	3mL Sterile Syringe Tray	Luer-Lock	25	8	200/CS
8881516200	6mL Sterile Syringe Tray	Luer-Lock	25	8	200/CS
8881512258	12mL Sterile Syringe Tray	Luer-Lock	25	8	200/CS
8881520251	20mL Sterile Syringe Tray	Luer-Lock	25	8	200/CS
8881535101	35mL Sterile Syringe Tray	Luer-Lock	12	10	120/CS
8881560232	60mL Sterile Syringe Tray	Luer-Lock	12	10	120/CS

BD Single-Use Syringe Trays



- Help improve workflow by reducing the number of single-syringe packaging that needs to be opened
- Support compliance with USP 797 packaging wipe-down requirement

Order No.	Description	Tip Type	Syringes per Tray	Trays per Case	Qty
309628	1mL Sterile Syringe Tray	Luer-Lock	25	12	300/CS
309702	3mL Sterile Syringe Tray	Luer-Lock	25	12	300/CS
309703	5mL Sterile Syringe Tray	Luer-Lock	25	12	300/CS
309605	10mL Sterile Syringe Tray	Luer-Lock	20	15	300/CS
305617	20mL Sterile Syringe Tray	Luer-Lock	10	20	120/CS
305618	30mL Sterile Syringe Tray	Luer-Lock	10	20	120/CS
309680	50mL Sterile Syringe Tray	Luer-Lock	20	10	120/CS
8881135609	35mL Plungerless Syringe (Syringe Only - No Tray)	Luer Lock	10	15	150/CS

ISO Certified Pharmaceutical & Bio-Pharma Filters

Talk to your sales rep for additional Inlet/Outlet fitting availability.
(Ex. Luer, Hose Barb, and Tri-Clamp Fittings)



Sterile Compounding PES Disk Filters



Order No.	Size	Layers	Qty
D25CS020LFLM-PH-ETO	25mm	Single	1/EA
D40CS020LFLM-PH-ETO	40mm	Single	1/EA
D65RS020LFLM-PH-ETO	65mm	Single	1/EA
D65RDS065S020LFLM-PH-ETO	65mm	Dual	1/EA
D90RS020LFLM-PH-ETO	90mm	Single	1/EA
D90RDS065S020LFLM-PH-ETO	90mm	Dual	1/EA

JKP Series - Small Capsule PES Filters



Order No.	Layers	Qty
JKPS020LFLM-ETO	Single	1/EA
JKPDS065S020LFLM-PH-ETO	Dual	1/EA

MKP Series - Small Capsule PES Filters



Order No.	Layers	Qty
MKPS020DLFLM-ETO	Single	1/EA
MKPDS065S020DLFLM-ETO	Dual	1/EA

SKL-H Series - Medium Capsule PES Filters with Female Luer Lock (Inlet) / Male Luer Lock (Outlet)



Order No.	Layers	Qty
SKLS020HLFLM-PH-ETO	Single	1/EA
SKLDS065S020HLFLM-PH-ETO	Dual	1/EA

SKL-S Series - Medium Capsule PES Filters with Female Luer Lock (Inlet) / Male Luer Lock (Outlet)



Order No.	Layers	Qty
SKLS020SLFLM-PH-ETO	Single	1/EA
SKLDS065S020SLFLM-PH-ETO	Dual	1/EA

SKL-L Series - Medium Capsule PES Filters with Hose Barb Fittings



Order No.	Layers	Qty
SKLDS065S020L2H2H-PH-ETO	Dual	1/EA



HCT PPE is White Low-Linting Fabric, Sterile, & Individually Packaged

Range:	Tolerance:
0-30cm	+/- 0.635cm
30-61cm	+/- 1.27cm
>61cm	+/- 2.54cm


Hoods

White 50GSM SF fabric. Features a tie-back to properly fit any user's head size.

Item ID:

- HCTHDIR

Size: One Size Fits All



Back of Hood

Coveralls


White 55GSM SF fabric.

Item ID:

- HCTCVRLIR-S
- HCTCVRLIR-M
- HCTCVRLIR-L
- HCTCVRLIR-XL
- HCTCVRLIR-2XL
- HCTCVRLIR-3XL
- HCTCVRLIR-4XL
- HCTCVRLIR-5XL

Size:

- Small
- Medium
- Large
- XL
- 2XL
- 3XL
- 4XL
- 5XL



Coveralls have zipper closure

Tunnel-knit cuffs on wrist with thumb loop

Boot Covers

White 63GSM SF fabric boots with anti-skid bottoms.

Item ID:

- HCTBOOTCIR-S/M
- HCTBOOTCIR-L/XL
- HCTBOOTCIR-2XL/3XL

Size:


- Small/Medium
- L/XL
- 2XL/3XL

Tunneled elastic on top

Double ties

Tunneled elastic on ankles

Back of coveralls



Tunneled elastic on back half of waist

HCT PPE Sizing Charts



Hoods (One Size Fits All)

Item ID: **HCTHDIR** | Qty: **100/CS**


SGS HCTCVRLIR Coverall Testing Results:

SGS Test Report #SL52405251581601TX

Composition: PP 60% + PE 40%

Testing Period: Jan 22, 2024 - Jan 31, 2024

Requirements Met: **YES**




Coveralls

Item ID: **HCTCVRLIR** | Qty: **25/CS**

Sizing	S	M	L	XL	2XL	3XL	4XL	5XL
Overall length	150	156	162	168	174	180	186	202
Torso length	88	90	92	94	98	102	106	112
Zipper length	70	70	70	75	75	78	78	78
Inseam	62	66	70	74	76	78	80	80
Chest width	92	100	108	116	124	132	140	148
Collar seam to waist elastic	47	48	54	55	56	57	59	61
Overall arm length	162	167	172	180	185	190	195	200
Sleeve length	55	56	58	60	60	64	64	66
Sleeve width	33	33	33	34	34	36	38	38
Relaxed wrist opening	8	8	8	8	8	8	8	8
Stretched wrist opening	15.2	15.2	16	16	16	16	17	17
Relaxed ankle opening	12	12	12	12	14	14	14	14
Stretched ankle opening	24.1	24.1	26	26.7	26.7	26.7	28	28
Relaxed waist opening	29	29	35	35	38	38	38	38
Stretched waist opening	52	52	60	60	62	62	62	64

Boot Covers

Item ID: **HCTBOOTCIR** | Qty: **50 Pairs/CS**



Sizing	S/M	L/XL	2XL/3XL	
Foot opening width	A	28	30	32
Boot height	B	41	46	51
Boot width	C	35	37	39
Anti-skid bottom width	D	13	13	14
Anti-skid bottom length	E	29	34	37
Upper tie length	F	29	36	45
Lower tie length	G	42	52	65
Tie width	H	2	2	2
Upper tie - Foot opening	I	5	5	5
Lower tie - Heel	J	10	10	10

Measurements are in centimeters

Essential Supplies for Your Pharmacy Environment

Face Protection

Order No.	Description	Qty
HCTWB24	Bouffant - 24" White	500/CS



PY.GB9410STM	Pyramex Latex Free Wrap-Around Goggles Dual Frame Design with Synthetic Rubber Gasket Anti-Fog, Protects from Splash and Dust H2MAX (Apply goggles onto face prior to disinfectant mask)	1/EA
--------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------



MO.4820	Moldex 4820 N95 Plus Relief Face Mask SmartStrap® & Full Foam Face Seal, Latex Free Metal-free, 100% PVC-Free	8/CS
---------	---------------------------------------------------------------------------------------------------------------------	------



Cardinal Heath™ ESTEEM™ CP Sterile Polyisoprene Gloves

- 1 pair per inner high-density polyethylene (PE) wallete
- Cuff Style: Beaded, interlocking
- Typical length: 12 in. / 305 mm
- Sterilization: Sterile – Gamma irradiated



Order No.	Size	Pouches per Bag	Outer Bags per Case	Qty
2Y72PL55	5½	50	4	200/CS
2Y72PL60	6	50	4	200/CS
2Y72PL65	6½	50	4	200/CS
2Y72PL70	7	50	4	200/CS
2Y72PL75	7½	50	4	200/CS
2Y72PL80	8	50	4	200/CS
2Y72PL85	8½	50	4	200/CS
2Y72PL90	9	50	4	200/CS



QI Medical Surface/Fingertip Testing

QI Medical EnviroTest™



- Ideal for use in critical surface and gloved fingertip sampling to meet USP <797> requirements. Unique pedestal cap and hinged paddle
- Stores at room temperature with long shelf life without refrigeration
- Log sheet & gummed labels included




EnviroTest™ TSA (trypticase soy agar) Growth Media Paddles

Order No.	Description	Shelf Life/Months	Qty
 ET1000	Double sided TSA (trypticase soy agar) growth media paddles for surface, air, or glove fingertip sampling, formulated to inactivate many antimicrobial compounds	12	10 Units/CS
 ET3000	Double sided MEA (malt extract agar) growth media paddles for surface, air, or glove fingertip sampling, formulated to inactivate many antimicrobial compounds with added antibiotics to prevent bacterial overgrowth	12	10 Units/CS

GroMed® Plates - Formulated with USP required TSA with Poly 80 & Lecithin

Order No.	Description	Shelf Life/Months	Qty
 GM0601	60mm raised (convex) surface contact plate (Shelf stable between 2-25°C)	4	10/Sleeve
 GM0901	90mm dish (concave) settling plate (Shelf stable between 2-25°C)	6	10/Sleeve

GroMed™ Incubators - Lab Cell Culture Incubator for Pharmacy or Lab Use

Order No.	Description	Qty
 QL10-140	Incubator 13"Wx11"Dx15"H, 3 shelves, dial thermometer	1/CS
QL12-140	Large incubator 19"Wx13"Dx21"H, 3 shelves, dial thermometer	1/CS

*Shelf life is based on the date of manufacture.

GroMed™ Growth Media, Test Kits, Vials, & More



Order No.	Description	Shelf Life/Months	Qty
 GM0030	TSB Growth media, 3mL of standard strength in glass ampules	24	20 Ampules/CS
 GM0100	TSB Growth media, 10mL tube, screw cap and needle access port	24	20 Tubes/CS
 GM0101	FTM Growth media, 10mL tube, screw cap and needle access port	12	20 Tubes/CS
 GM0200	TSB Growth media, 20mL vial, needle/spike septum	24	20 Vials/CS
 GM0201	FTM Growth media, 20mL vial, needle/spike septum	12	10 Vials/CS
 GM0203	Fluid D, 20mL vial, needle/spike septum	12	20 Vials/CS
 GM0252	TSB Growth media, 25mL vial, 2X, needle septum	24	10 Vials/CS
 GM1400	TSB Growth media, 100mL bag with needle access port and spike port	24	10 Bags/CS
 GM1420	TSB Growth media, 100mL bag with needle access port and spike port	24	10 Bags/CS
 GM1500	TSB Growth media, 100mL bag with Male Luer Lock port and spike port	24	10 Bags/CS
 GM2000	TSB Growth media, 100mL vial with septum for needle or dispensing pin access	24	10 Vials/CS

*Shelf life is based on the date of manufacture.

QI Medical GroMed™ Aseptic Technique Training Aids



Order No.	Description	Shelf Life/Months	Qty
 GM2001	FTM Growth media, 100mL vial with septum for needle or dispensing pin access	24	10 Vials/CS
 GM3000	Soybean-Casein Digest Media, approximately 3 grams of nonsterile powder in foil pouch, makes 100mL of TSB media	48	20 Pouches/CS
 GM4000	PDA Potato dextrose agar, 10mL tubes, test for fungal contamination	24	20 Tubes/CS
 GM5000	TSB Growth media, 500mL bag, 1X strength, breakaway covered needle access port and spike port	24	10 Bags/CS
 GM5500	TSB Growth media, 1L bag with needle access port and spike port	24	5 Bags/CS
 GM5900	TSB Growth media, 4L bag with needle access port and spike port	24	2 Bags/CS
 GM6005	5cc TSB Growth media in 12cc luer, lock syringe	12	10 Syringes/CS
 GM6105	5cc FTM (Fluid Thioglycollate) Growth Media in 12cc luer lock syringe	24	10 Syringes/CS
 GM7020	PATT™, Personal Aseptic Technique test kit, 20mL vials, 100mL partially filled minibags, test procedure and results log	24	5/Vials - 5 Bags
 GM7030	PATT2™, Personal Aseptic Technique test kit, 3mL ampules, 20mL vials, 100mL partially filled minibags, test procedure and results log	24	5 Ampules 5 Vials - 5 Bags
 EV0200	20mL sterile empty vial for use in aseptic technique tests. Not for patient use	36	20 Vials/CS
 EV2000	100mL sterile empty vial for use in aseptic technique tests. Not for patient use	36	10 Vials/CS

*Shelf life is based on the date of manufacture.






QI Medical Contamination Testers



TuffTest2™ System - Aliquot Test for Emulsions and Suspensions

Order No.	Description	Shelf Life/Months	Qty
 GM8100	TuffTest2™ contamination tester, 100mL TSB with TSA in sterile plastic bottle, needle access port	12	4 Testers/CS
 GM8101	TuffTest2™ contamination tester, 100ml FTM with Clostridial agar in sterile plastic bottle, needle access port	12	4 Testers/CS

QuickTest™ System - Full Filtration for All LVPs, Including Emulsions

Order No.	Description	Shelf Life/Months	Qty
 QT1000	Contamination tester for use with: <ul style="list-style-type: none"> • ADComp™ Compounder • Automix Compounder (Clintec) • Autocomp 6™ Compounder • Exacta-Mix™ Compounder (Baxa) • Nutrimix™ Compounder (Abbott) 	24	10 Testers 10 100mL Bags of media
 QT1100	Contamination tester with universal spike and 18" extension set	24	10 Testers 10 100mL Bags of media 10 16 ga. Needles
 QT2000	Contamination tester for use with: <ul style="list-style-type: none"> • Any pump set with male luer connector • Repeater™ Pump (Baxa) • Hyperformer™ Compounder (McGaw) • PharmaAssist™ Pump (Excelsior) 	24	10 Testers 10 100mL Bags of media
 QT3000	Contamination tester for use with: Pinnacle compounder (BBraun)	24	10 Testers 10 100mL Bags of media
 QT4000	Contamination Tester for use with: APEX Compounder	24	10 Testers 10 100mL Bags of media


*Shelf life is based on the date of manufacture.

QI Medical Testing Systems & More



QT Micro™ Systems - Full Filtration for Small Volume Sterile Solutions

Order No.	Description	Shelf Life/Months	Qty
 TM6000	0.22µ Contamination tester, female luer lock inlet, female luer lock outlet, TSB Media	24	10 Testers 10 5cc Syringes of media
 TM6001	0.22µ Contamination tester, female luer lock inlet, female luer lock outlet, FTM Media	12	10 Testers 10 5cc Syringes of media
 TM6100	0.22µ Contamination tester, female luer lock inlet, male luer slip outlet, TSB Media	24	10 Testers 10 5cc Syringes of media
 TM6101	0.22µ Contamination tester, female luer lock inlet, male luer slip outlet, FTM Media	12	10 Testers 10 5cc Syringes of media
 GM0203	Fluid D, 20mL vial, needle/spike septum	12	20 Vials

ChemoTest™ - Hazardous Drug Handling Test Kits

Order No.	Description	Shelf Life/Months	Qty
 CT9530	Validate the technique of health care workers who manipulate hazardous drugs. Uses fluorescein dye and a UV light to detect breaks in safe technique and effectiveness of hazmat clean up procedures. 3mL and 50mL vials	24	5 Tests

PALL Syringe Filters















Order No.	Description	Qty
 HP1002	0.22µ filter, 25mm, sterile, female luer lock inlet, male luer lock outlet, hydrophilic	50 Filters/CS
 PG005	Bubble Point Tester for use with syringe filters, 0-100psi range	1 Tester/CS

*Shelf life is based on the date of manufacture.



Contec Wipes Designed that Exceed USP Standards















Order No.	Description	Size	Qty
 HCPS6044	Contec Critical Site® Sterile wipes	4"x4"	1440 wipes/CS (60/pouch)
 HCPS2328IR	PROSAT® Sterile presaturated meltblown polypropylene wipes with 70% IPA/30% DI water	9"x11"	810 wipes/CS (30/pouch)
 HCPS0002IR	PROSAT® Sterile presaturated cellulose/polyester wipes with 70% IPA/30% DI water	9"x11"	810 wipes/CS (30/pouch)
 HCPS2328	PROSAT® presaturated meltblown polypropylene wipes with 70% IPA/30% DI water	9"x11"	810 wipes/CS (30/pouch)
 HCPS0002	PROSAT® presaturated cellulose/polyester wipes with 70% IPA/30% DI water	9"x11"	810 wipes/CS (30/pouch)
 PS-911EB	PROSAT® Sterile presaturated meltblown polypropylene wipes with 70% IPA/30% DI water	9"x11"	1440 wipes/CS (30/pouch)
 PS-911	PROSAT® presaturated meltblown polypropylene wipes with 70% IPA/30% DI water	9"x11"	1440 wipes/CS (30/pouch)
 PS-911-12	PROSAT® presaturated meltblown polypropylene wipes with 70% IPA/30% DI water	9"x11"	360 wipes/CS (30/pouch)
 PSPS0076	PROSAT® Sterile Polynit Heatseal Half-Fold Wipes with 70% IPA/30% DI water	9"x9"	800 wipes/CS (20/pouch)
 HCTB3315IR	Contec Healthcare TB1-3300™ disinfectant Sterile wipes	9"x11"	810 wipes/CS (30/pouch)
 HCTB3315	Contec Healthcare TB1-3300™ disinfectant wipes	9"x11"	810 wipes/CS (30/pouch)
 PSC20005	PROSAT® Sterile Theta™ wipes, lightweight nonwoven cellulose/polyester presaturated 70% IPA/30% DI water	9"x11"	1400 wipes/CS (28/pouch)
 21221	Contec PREempt RTU disinfectant, meltblown polypropylene wipes with Accelerated Hydrogen Peroxide	6"x7"	12/CS (160/canister)
 2B221	PREempt Plus Meltdown Wipes, polypropylene wipes with Accelerated Hydrogen Peroxide	6"x7"	12/CS (160/canister)











Contec Wipes (cont.)




Contec Low-Lint, Sterile/Non-Sterile Dry Wipes



Order No.	Description	Size	Qty
 HCC20001	Absorbent, disposable, low-linting Hand Drying wipes.	17" x 9"	900 wipes/CS (60/pouch)
 9451	Wall mounted stainless steel dispenser for HCC20001 hand drying wipes	N/A	1/CS
 NWPZ0001	Amplitude™ Kappa™ Sterile LE wipes lyocell/polyester blend nonwoven fabric	9" x 9"	1000 wipes/CS (30/pouch)
 NWPZ0002	Amplitude™ Kappa™ Sterile LE wipes lyocell/polyester blend nonwoven fabric	12" x 12"	1000 wipes/CS (30/pouch)
 AMEP0001	Amplitude™ Epsilon™ wipes cellulose/polyester nonwoven wipes	9" x 9"	3600 wipes/CS (30/pouch)
 AMEP0002	Amplitude™ Epsilon™ wipes cellulose/polyester nonwoven wipes	12" x 12"	3600 wipes/CS (30/pouch)
 AMEC0002	Amplitude™ EcoCloth™ wipes cellulose/polyester nonwoven wipes	6" x 6"	6000 wipes/CS (30/pouch)
 AMEC0003	Amplitude™ EcoCloth™ wipes cellulose/polyester nonwoven wipes	9" x 9"	3600 wipes/CS (30/pouch)
 AMEC0004	Amplitude™ EcoCloth™ wipes cellulose/polyester nonwoven wipes	12" x 12"	2100 wipes/CS (30/pouch)
 AMEC0005	Amplitude™ EcoCloth™ wipes cellulose/polyester nonwoven wipes	18" x 18"	750 wipes/CS (10/pouch)
 AMSI0001	Amplitude™ Sigma™ Wipes Cellulose/polyester hydroentangled nonwoven wipes	9" x 9"	3600 wipes/CS (30/pouch)
 C2-99IR/25	SterileSorb™ wipes Sterile hydroentangled cellulose/polyester wipes	9" x 9"	1050 wipes/CS (30/pouch)
 C2-1212IR	SterileSorb™ wipes Sterile hydroentangled cellulose/polyester wipes	12" x 12"	1050 wipes/CS (30/pouch)

Order No.	Description	Size	Qty
 LWPS0006	Polynit Heatseal Sterile wipes textured polyester interlock knit	9"x9"	1200 wipes/CS (25/pouch)
 LWMM0300	Microfiber PET/PP wipes Cleanroom laundered, bulk	10"x18"	400 wipes/CS (100/pouch)
 LWMM0400	Microfiber PET/PP wipes Cleanroom laundered, bulk	10"x18"	1 kit/CS
 HCPW3001	Microfiber PET/PP wipes Cleanroom laundered, bulk	10"x12"	700 wipes/CS (25/pouch)
 HCPW3003	Microfiber PET/PP wipes Cleanroom laundered, bulk	10"x18"	250 wipes/CS (25/pouch)
 85337	Sanotex® environmental surface wipes Half fold	10"x12"	1800 wipes/CS (60/roll)
 85338	Sanotex® environmental surface wipes Half fold	10"x12"	320 wipes/CS (40/roll)
 85339	Sanotex® environmental surface wipes	6"x8"	900 wipes/CS (150/roll)
 85340	Sanotex® environmental surface wipes Half fold	10"x12"	700 wipes/CS (10/roll)
 85341	Sanotex® environmental surface canister Empty self-closing dispenser, reusable	N/A	6 canisters/CS



Contec Critical Pharmacy Products







Contec Cleaning and Disinfection Chemicals

Specially-formulated sporicidal agents, disinfectants, and chemicals for cleaning and disinfecting sterile/critical pharmacy environments; including 70% Isopropanol, hydrogen peroxide-based solutions, to fast-acting sporicidal disinfectants and broad-spectrum daily disinfectants.



Order No.	Description	Size	Qty
 HC85335	PeridoxRTU broad spectrum sporicidal disinfectant EPA Reg. No. 8383-13 Flip-top cap	32 oz.	6/CS
 HC85335IR	PeridoxRTU, Sterile broad spectrum sporicidal disinfectant EPA Reg. No. 8383-13, Flip-top cap	32 oz.	6/CS
 HC85335IR	PeridoxRTU broad spectrum sporicidal disinfectant EPA Reg. No. 8383-13	1 gal.	6/CS
 HC85336IR	PeridoxRTU, Sterile broad spectrum sporicidal disinfectant EPA Reg. No. 8383-13	1 gal.	6/CS
 HCTB3315IR	Contec Healthcare TB1-3300 disinfectant Sterile wipes	9"x11"	810 wipes/CS (30/pouch)
 HCTB3315	Contec Healthcare TB1-3300 disinfectant wipes	9"x11"	810 wipes/CS (30/pouch)
 HCTB3300	Contec Healthcare TB1-3300 disinfectant Single bagged, push-pull cap	32 oz.	12/CS
 HCTB3300IR	Contec Healthcare TB1-3300 disinfectant Sterile , double bagged, push-pull cap	32 oz.	12/CS
 HCTB3305	Contec Healthcare TB1-3300 disinfectant Single bagged	1 gal.	4/CS
 HCTB3305IR	Contec Healthcare TB1-3300 disinfectant Sterile , double bagged	1 gal.	4/CS
 SB167030IR	Contec Sterile 70% isopropanol Trigger sprayer	16 oz.	12/CS
 SB327030IR	Contec Sterile 70% isopropanol Trigger sprayer	32 oz.	12/CS
 SB1287030IR	Contec Sterile 70% isopropanol	1 gal.	4/CS

Contec Cleaning and Disinfection Chemicals (cont.)






Order No.	Description	Size	Qty
 HCFT7030IR	Contec Sterile 70% isopropanol Flip-top cap	16 oz.	12/CS
 HCFT7030IR-32	Contec Sterile 70% isopropanol Flip-top cap	32 oz.	12/CS
 21101	Contec PREempt RTU Disinfectant With Accelerated Hydrogen Peroxide, capped	32 oz.	12/CS
 21105	Contec PREempt RTU Disinfectant With Accelerated Hydrogen Peroxide, capped	1 Gallon	4/CS
 2B101	Contec PREempt Plus One-Step Disinfectant Cleaner	32 oz.	12/CS
 2B105	Contec PREempt Plus One-Step Disinfectant Cleaner	1 Gallon	4/CS

Contec Low-Lint Mops for Cleanrooms & More

Systems and accessories for low-lint mopping that meet and exceed USP Standards, including mop frames, handles, wall washing systems, single-use disposable microfiber mop pads, and more









Contec Cleaning Pads / Mop Covers (One-Time Use, Disposable)















Order No.	Description	Size	Qty
 PSME0001	EasyReach™ cleaning pads Quilted polyester presaturated with 70% IPA	2.75" x 7.75"	64/CS
 HCMT0011	EasyReach™ cleaning pads with ultrasonic pockets. Quilted polyester	2.75" x 7.75"	60/CS
 HCMT0015	EasyReach™ cleaning pads with ultrasonic pockets. Quilted polyester, Sterile	2.75" x 7.75"	60/CS
 HCPM5001	MicroCinch™ Mop Microfiber mop cover with pockets. Standard weight	5" x 15.5"	120/CS
 HCPM5051	MicroCinch™ Mop Microfiber mop cover with pockets. Standard weight	5" x 15.5"	70/CS

Contec Mops & Handles

Contec Mop Heads (One-Time Use, Disposable)

Order No.	Description	Size	Qty
 HCKM3001	Klean Max™ mop head, small	8"x4"x1"	48/CS
HCKM3002	Klean Max™ mop head, medium	13"x4"x1"	28/CS
HCKM3003	Klean Max™ mop head, large	15"x4"x1"	24/CS
HCKM3051	Klean Max™ mop head, small. Sterile	8"x4"x1"	48/CS
HCKM3052	Klean Max™ mop head, medium. Sterile	13"x4"x1"	28/CS
HCKM3053	Klean Max™ mop head, large. Sterile	15"x4"x1"	24/CS
 HCKM3005	Klean Max™ mop head, small. Sealed edge	9.5"x5"x0.375"	40/CS
HCKM3004	Klean Max™ mop head, medium. Sealed edge	15"x5"x0.375"	16/CS
HCKM3055	Klean Max™ mop head, small Sealed edge, Sterile	9.5"x5"x0.375"	40/CS
HCKM3054	Klean Max™ mop head, medium Sealed edge, Sterile	15"x5"x0.375"	24/CS
 EZMP0350	EasyCurve™ mop head. Microfiber individually packaged, Sterile	14"x4.3"x.75"	1/CS
EZMP0300	EasyCurve™ mop head. Microfiber	14"x4.3"x.75"	1/CS
 TASK0100	Quicklean foam mop head Laminated with polyester fabric	N/A	32/CS
 HCPM2002	Premira II White Microfiber Pads	N/A	200/CS
 BLVKMOPIR	EasySat™ Mop Head, Sterile Knitted polyester/foam mop	21.5"x18.25"x4.75"	28/CS

Contec Mop Handles / Mop Frames / Tools

Order No.	Description	Size	Qty
 2642	QuickConnect™ anodized aluminum handle (Compatible with EasyReach™)	16"-30"	1/CS
 2643	QuickConnect™ anodized aluminum handle (Compatible with EasyReach™)	50"-92"	1/CS
 2644	QuickConnect™ stainless steel handle (Compatible with EasyReach™ & EasyCurve™)	16"-30"	1/CS
 2645	QuickConnect™ stainless steel handle (Compatible with EasyCurve™)	50"-92"	1/CS
 2725E	QuickConnect™ stainless steel handle (Compatible with MicroCinch™, EasyCurve™, Klean Max™, & QuickKlean™)	60"	1/CS
 2746S	QuickConnect™ handle Fiberglass with stainless steel tip	40"-70"	1/CS
 2685LF	QuickConnect™ EasyReach™ mop frame, layflat. Stainless steel	7.5"	1/CS
 2665LF	EasyReach™ Cleaning Tool, layflat, lightweight Stainless steel	16.5"	1/CS
 HCPH8007	MicroCinch™ flat hinged mop frame Stainless steel	16"x4.25"	1/CS
 2651SS	QuickConnect™ mop frame, small For Klean Max™ small mop head	7.5"	1/CS
 2652LF	QuickConnect™ mop frame, layflat, medium For Klean Max™ medium and large mop heads	11.5"	1/CS
 2740	MicroCinch™ flat hinged mop frame Blue polypropylene	16"x4.25"	1/CS
 2745	QuickKlean flat mop head frame Blue polypropylene	16"	1/CS
 2724	EasyCurve™ pivoting mop frame Stainless steel	2.75"x5.25"x14.5"	1/CS

Contec Garb & PPE



Cleanroom garb for critical and USP-compliant environments, including gowns, frocks, garb, sterile personal protective equipment (PPE), and more



Contec CritiGear™ Sterile Nitrile Gloves (One-Time Use, Disposable)

- Contec CritiGear Sterile Nitrile Gloves have been specifically designed for pharmacy environments regulated by USP <797> and USP <800>.
- Sterility Assurance Level 10^{-6} | Tested per ASTM D6978-05

Order No.	Description	Packaging	Qty
HCGS0002	X-Small	1 pair/pack	200 pair/CS
HCGS0012	Small	1 pair/pack	200 pair/CS
HCGS0022	Medium	1 pair/pack	200 pair/CS
HCGS0032	Large	1 pair/pack	200 pair/CS
HCGS0042	X-Large	1 pair/pack	200 pair/CS
HCGS0052	2X-Large	1 pair/pack	200 pair/CS

Contec Shoe Covers (One-Time Use, Disposable)

Order No.	Description	Size	Qty
 SCR645L-150	SCR645 series shoe covers. Blue	Universal	300 each/CS
SCR645XL-150	SCR645 series shoe covers. Blue	X-Large	300 each/CS
SCR645XXL-150	SCR645 series shoe covers. Blue	2X-Large	300 each/CS
 SCPL500	SCPL series compressed polyethylene shoe covers. White	Universal	1000 each/CS
 LKBT0001-S/M-IR	CleanMax™ STERILE 19" Boot Cover Non-Skid PVC vinyl sole, tunneled elastic top	Small/Medium	50/CS



Bouffant Cap 24" - White

- Strong, Durable, & Lightweight with No Latex Proteins
- Contamination Control
- Individually Inspected for Stitching, Sizing, Material and Workmanship

Order No.	Description	Qty
NBC2400	24"	1000/CS



Beard Cover White Spunbonded Poly - White

- Strong, Durable, & Lightweight with No Latex Proteins
- Contamination Control
- Individually Inspected for Stitching, Sizing, Material and Workmanship

Order No.	Description	Qty
HBC1000	Approx. Length: 38.1cm (15 in.)	1000/CS



Bouffant Cap 24" - White

- Strong, Durable, & Lightweight with No Latex Proteins
- Contamination Control
- Individually Inspected for Stitching, Sizing, Material and Workmanship

Order No.	Description	Qty
GOGLIR	24"	1000/CS



Contec Cleanroom Mask - White

- Earloop Face Mask
- Strong, Durable, & Lightweight
- Cleanroom Rated

Order No.	Description	Qty
MAS-9010	White	500/CS

Contec Garb & PPE



Contec CritiGear™ Cleanroom Frock with Thumb Loop

Order No.	Size	Qty
HCGA0002	X-Small	30/CS
HCGA0012	Small	30/CS
HCGA0022	Medium	30/CS
HCGA0032	Large	30/CS
HCGA0042	X-Large	30/CS
HCGA0052	2X-Large	30/CS
HCGA0062	3X-Large	30/CS
HCGA0072	4X-Large	30/CS
HCGA0082	5X-Large	30/CS



Contec CritiGear™ Chemo Hazardous Drug Gowns with Thumb Loop

Order No.	Size	Qty
HCGA0800-S/M	Small/Medium	30/CS
HCGA0800-L/XL	Large/X-Large	30/CS
HCGA0800-2XL/3XL	2X-Large/3X-Large	30/CS
HCGA0800-4XL/5XL	4X-Large/5X-Large	30/CS
HCGA0800-6XL/7XL	6X-Large/7X-Large	30/CS



CleanMax™ Coveralls with Hood, Tunneled Wrists, & Ankles

Order No.	Size	Qty
CVHD0022	Medium	25/CS
CVHD0032	Large	25/CS
CVHD0042	X-Large	25/CS
CVHD0052	2X-Large	25/CS
CVHD0062	3X-Large	25/CS

STERILE



ClearMax™ garments feature bound seams, which are sewn with an additional outer binding to increase seam strength and provide better barrier protection.

CleanMax™STERILE Coveralls with Hood, Tunneled Wrists, & Ankles

Order No.	Size	Qty
CVHD0022IR	Medium	25/CS
CVHD0032IR	Large	25/CS
CVHD0042IR	X-Large	25/CS
CVHD0052IR	2X-Large	25/CS
CVHD0062IR	3X-Large	25/CS



CVRL

CleanMax™ Coveralls with Tunneled Wrists & Ankles

Order No.	Size	Qty
CVRL0022	Medium	25/CS
CVRL0032	Large	25/CS
CVRL0042	X-Large	25/CS
CVRL0052	2X-Large	25/CS
CVRL0062	3X-Large	25/CS

STERILE



CVRL

CleanMax™STERILE Coveralls with Tunneled Wrists & Ankles

Order No.	Size	Qty
CVRL0022IR	Medium	25/CS
CVRL0032IR	Large	25/CS
CVRL0042IR	X-Large	25/CS
CVRL0052IR	2X-Large	25/CS
CVRL0062IR	3X-Large	25/CS



Contec Critical Environment Supplies

Critical supplies for your critical environment, including prep and tack mats, disposable mesh foam sponges, and an essential hazardous drug spill kit



Clean Scrub™ Sponge (Disposable mesh sponge with foam core)

- Textured polyester mesh removes heavily soiled contamination from surfaces
- Compatible with a wide range of disinfectants

Order No.	Size	Qty
VKSC0004	5" x 3" x 0.6"	48 each/CS



Contec Healthcare Essential HD Spill Toolkit (for up to 1000mL)

- Activated charcoal spill pads begin chemical deactivation of spills, absorbs, and renders liquid drug unrecoverable, enhancing the safety of workers absorbing up to 1000mL per pad
- Integrated Chemo Bio-Wipe bag safely contains macro spill components for later disposal in RCRA* container
- Plastic scoop and rigid container conveniently facilitate powder spill cleanup and containment

Order No.	Qty
HCSK0001	2 kits/CS



Irradiated Bucket Liners

- 4 mil polyethylene; fits 5 gallon, 6.5 gallon, 7 gallon
- 2/bag, 20 bags/case


Order No.	Size	Qty
BKTBAG	34" x 24" x .004"	40/CS



Contec Prep Mat - Contaminaton Control Mat

- Low-lint construction; absorbent and barrier layers
- Lot number on inner bags for traceability
- Available irradiated and nonsterile

Order No.	Description	Size	Qty
HCPD0001	Contec prep mat	16"x21"	600 each/CS
HCPD0005	Contec prep mat. Irradiated	16"x21"	600 each/CS



Contec Tack Mat - Contaminaton Control Flooring

- Placed immediately before the entrance to controlled areas to reduce the introduction of contamination
- Thirty peelable layers for convenience
- Multiple sizes are available
- Frames are an optional way to avoid exposing the floor to adhesives

Order No.	Description	Size	Qty
TMF2436	Tacky mat frame	24"x36"	1 each/CS
TMF1836	Tacky mat frame	18"x36"	1 each/CS
TM0001	Tacky mat. Blue	24"x36"	120/CS
TM0006	Tacky mat. Blue	18"x36"	120/CS



Medipak UVLI® Covers and Bags for UV Light-Sensitive Drug Protection



- Using Genuine UVLI® Covers and Bags is a cost-effective solution to protect liquid drugs administered intravenously, or solid drug transit from the harmful effects of UV light
- Available in 13 different sizes
- Allow for the administrator to read through the bag, eliminating the need to lift the cover to check levels and labels
- Talk to your sales rep about product programs and availability




Available Colors




Semi-transparent **Amber** UVLI™ covers make it easy to read labels and levels



Semi-transparent **Green** UVLI™ covers make it easy to read labels and levels



Black UVLI™ covers block 98% or more of light transmission and are nearly opaque



Brown UVLI™ covers block 99% or more of light transmission and are opaque

UVLI Tubing Covers for IV Lines


Order No.	Color	Size	Qty
0590	Amber	2.5" x 2500'	1 Roll/CS
590X	Amber	3" x 72"	100 Precut/CS
0591	Amber	2.5" x 1000'	1 Roll/CS
0593	Amber	2.5" x 300'	1 Roll/CS
0790	Green	2.5" x 2500'	1 Roll/CS
790X	Green	3" x 72"	100 Precut/CS
0791	Green	2.5" x 1000'	1 Roll/CS
0793	Green	2.5" x 300'	1 Roll/CS
0890	Black	2.5" x 2500'	1 Roll/CS
870X	Black	3" x 72"	100 Precut/CS
0891	Black	2.5" x 1000'	1 Roll/CS
0893	Black	2.5" x 300'	1 Roll/CS
0990	Brown	2.5" x 2500'	1 Roll/CS
970X	Brown	3" x 72"	100 Precut/CS
0991	Brown	2.5" x 1000'	1 Roll/CS
0993	Brown	2.5" x 300'	1 Roll/CS



UVLI Zip Bags - Multi-Use Closure

UVLI™ Zip Bags provide a convenient multi-use closure, are cost-effective, and disposable for protecting light-sensitive medications during transportation and administration. Our bags exceed the United States Pharmacopoeia Standard blocking over 90% of harmful UV light between 290nm and 450nm, so you can protect with confidence.

Order No.	Color	Size	Qty
AZ35	Amber	3"x5"	1000 Bags/CS
GZ35	Green	3"x5"	1000 Bags/CS
AZ39	Amber	3"x9"	1000 Bags/CS
GZ39	Green	3"x9"	1000 Bags/CS
AZ46	Amber	4"x6"	1000 Bags/CS
GZ46	Green	4"x6"	1000 Bags/CS
AZ58	Amber	5"x8"	1000 Bags/CS
GZ58	Green	5"x5"	1000 Bags/CS
AZ69	Amber	6"x9"	1000 Bags/CS
GZ69	Green	6"x9"	1000 Bags/CS
AZ84	Amber	8"x14"	1000 Bags/CS
GZ84	Green	8"x14"	1000 Bags/CS
AZ92	Amber	9"x12"	1000 Bags/CS
GZ92	Green	9"x12"	1000 Bags/CS



Clear Poly Bags for Syringes (Medipak Clear Polybags do not protect against UV light)

Order No.	Color	Size	Qty
0210	Clear	2.5"x10.5"	1000/CS
0315	Clear	3"x15"	1000/CS
0610	Clear	2.5"x2"	1000/CS
0615	Clear	2.5"x5"	1000/CS
0620	Clear	2.5"x8.5"	1000/CS

Medipak UVLI® Covers

UVLI Regular Covers - Seal Completely Across One End

Amber Covers

Order No.	Intended Use	Size	Qty
0510	Convenient for transport of tablets	2.5"x2"	1000 Bags/CS
0515	For large ampules	2.5"x5"	1000 Bags/CS
0520	For syringes	2.5"x8.5"	1000 Bags/CS
0530	For large syringes	3"x14"	1000 Bags/CS
0540	For 450ml and smaller IV bags	5"x8.5"	500 Bags/CS
0550	For 1/2-liter IV bags (500 ml)	6"x10"	1000 Bags/CS
0555	For 1-liter IV bags (1000 ml)	6"x14"	1000 Bags/CS
0560	For 1-liter IV bags (1000 ml)	8"x14"	500 Bags/CS
0565	For 1-liter IV bags (1000 ml)	8"x14"	1000 Bags/CS
0580	For 3-liter IV bags (3000 ml)	10"x18"	250 Bags/CS
0585	For 4-liter IV bags (4000 ml)	12"x18"	250 Bags/CS

Black Covers

0820	For syringes	2.5"x8.5"	1000 Bags/CS
0830	For large syringes	3"x14"	1000 Bags/CS
0840	For 450ml and smaller IV bags	5"x8.5"	500 Bags/CS
0850	For 1/2-liter IV bags (500 ml)	6"x10"	1000 Bags/CS
0855	For 1-liter IV bags (1000 ml)	6"x14"	1000 Bags/CS
0860	For 1-liter IV bags (1000 ml)	8"x14"	500 Bags/CS
0865	For 1-liter IV bags (1000 ml)	8"x14"	1000 Bags/CS
0880	For 3-liter IV bags (3000 ml)	10"x18"	250 Bags/CS
0885	For 4-liter IV bags (4000 ml)	12"x18"	250 Bags/CS

Green Covers

Order No.	Intended Use	Size	Qty
0715	For large ampules	2.5"x5"	1000 Bags/CS
0720	For syringes	2.5"x8.5"	1000 Bags/CS
0730	For large syringes	3"x14"	1000 Bags/CS
0740	For 450ml and smaller IV bags	5"x8.5"	500 Bags/CS
0750	For 1/2-liter IV bags (500 ml)	6"x10"	1000 Bags/CS
0755	For 1-liter IV bags (1000 ml)	6"x14"	1000 Bags/CS
0760	For 1-liter IV bags (1000 ml)	8"x14"	500 Bags/CS
0765	For 1-liter IV bags (1000 ml)	8"x14"	1000 Bags/CS
0780	For 3-liter IV bags (3000 ml)	10"x18"	250 Bags/CS
0785	For 4-liter IV bags (4000 ml)	12"x18"	250 Bags/CS



Brown Covers

0920	For syringes	2.5"x8.5"	1000 Bags/CS
0930	For large syringes	3"x14"	1000 Bags/CS
0940	For 450ml and smaller IV bags	5"x8.5"	500 Bags/CS
0950	For 1/2-liter IV bags (500 ml)	6"x10"	1000 Bags/CS
0955	For 1-liter IV bags (1000 ml)	6"x14"	1000 Bags/CS
0960	For 1-liter IV bags (1000 ml)	8"x14"	500 Bags/CS
0965	For 1-liter IV bags (1000 ml)	8"x14"	1000 Bags/CS
0980	For 3-liter IV bags (3000 ml)	10"x18"	250 Bags/CS
0985	For 4-liter IV bags (4000 ml)	12"x18"	250 Bags/CS

UVLI Slit-Top Covers - Has Gap in the Seal for Easy IV Hook Use

Amber Covers

Order No.	Intended Use	Size	Qty
C581	For 3-liter IV bags (3000 ml)	10"x14"	250 Bags/CS
0541	For 450ml and smaller IV bags	5"x8.5"	500 Bags/CS
0551	For 1/2-liter IV bags (500 ml)	6"x10"	1000 Bags/CS
0556	For 1-liter IV bags (1000 ml)	6"x14"	1000 Bags/CS
0561	For 1-liter IV bags (1000 ml)	8"x14"	500 Bags/CS
0566	For 1-liter IV bags (1000 ml)	8"x14"	1000 Bags/CS
0581	For 3-liter IV bags (3000 ml)	10"x18"	250 Bags/CS
0586	For 4-liter IV bags (4000 ml)	12"x18"	250 Bags/CS

Black Covers

0841	For 450ml and smaller IV bags	5"x8.5"	500 Bags/CS
0851	For 1/2-liter IV bags (500 ml)	6"x10"	1000 Bags/CS
0856	For 1-liter IV bags (1000 ml)	6"x14"	1000 Bags/CS
0861	For 1-liter IV bags (1000 ml)	8"x14"	500 Bags/CS
0866	For 1-liter IV bags (1000 ml)	8"x14"	1000 Bags/CS
0881	For 3-liter IV bags (3000 ml)	10"x18"	250 Bags/CS
0886	For 4-liter IV bags (4000 ml)	12"x18"	250 Bags/CS

Green Covers

Order No.	Intended Use	Size	Qty
0741	For 450ml and smaller IV bags	5"x8.5"	500 Bags/CS
0751	For 1/2-liter IV bags (500 ml)	6"x10"	1000 Bags/CS
0756	For 1-liter IV bags (1000 ml)	6"x14"	1000 Bags/CS
0761	For 1-liter IV bags (1000 ml)	8"x14"	500 Bags/CS
0766	For 1-liter IV bags (1000 ml)	8"x14"	1000 Bags/CS
0781	For 3-liter IV bags (3000 ml)	10"x18"	250 Bags/CS
0786	For 4-liter IV bags (4000 ml)	12"x18"	250 Bags/CS

Brown Covers

0941	For 450ml and smaller IV bags	5"x8.5"	500 Bags/CS
0951	For 1/2-liter IV bags (500 ml)	6"x10"	1000 Bags/CS
0956	For 1-liter IV bags (1000 ml)	6"x14"	1000 Bags/CS
0961	For 1-liter IV bags (1000 ml)	8"x14"	500 Bags/CS
0966	For 1-liter IV bags (1000 ml)	8"x14"	1000 Bags/CS
0981	For 3-liter IV bags (3000 ml)	10"x18"	250 Bags/CS
0986	For 4-liter IV bags (4000 ml)	12"x18"	250 Bags/CS





Pharmacy Products from HCT



Orders will ship the same day if placed before 4 PM EST.
Talk to your sales rep about product programs and availability

To Learn More or Place an Order, Please Call Your Local Sales Rep:



Office: 508-881-6400
Fax: 508-881-3555
cs@hctsolutions.com
www.hctsolutions.com